



**VostraMed**  
Medical Technology

# NephroPatch™

External puncture closure

Selective compression to ensure functional efficiency of the shunt.

## Ordering Information:

| Reference Number (EAN) | Description              |
|------------------------|--------------------------|
| 7629999041318          | NephroPatch™ M SU (25x2) |

1 sales unit contains 25 pouches NephroPatch™

## NephroPatch™

NephroPatch™ is a trademark of VostraMed. Photo(s) on file at VostraMed. All drawings are artist's representations only and should not be considered as an engineering drawing or photograph. Procedures with above device shall only be carried out by physicians/medical professionals who are familiar with the use of the device as well as with possible complications. Prior to use, it is important to read package insert thoroughly for instructions for use, warnings and potential complications associated with the use of this device. Please check regulatory status of the device before distribution in areas where CE marking is not the regulation in force.

© 2017 VostraMed. All rights reserved.

10028-001/003-2019-10

### Distributor:

### For more Information:

VostraMed AG  
Baarer mattstrasse 6  
CH-6340 Baar  
Switzerland  
Phone: +41 41 552 25 10  
[www.vostramed.ch](http://www.vostramed.ch)

Pulsatile compression – quick – comfortable


1

*Perform dialysis in the usual manner. Leave fistula needle in situ. Expected skin adhesion sites have to be fat-free and dry.*



2

*Pull fistula needle 1 to 2 cm out of the shunt. Position pressure pad over the puncture channel. To perform as intended NephroPatch™ M must be attached circumferentially but holding strips must not overlap.*



3

*Promptly remove fistula needle and keep digital pressure for a little while.*



4

*NephroPatch™ M should be fastened in such a manner as to allow for haemostasis while continuous blood flow is maintained. Remove NephroPatch™ M after 4, preferably after 8–14 hours, in the same way as a plaster.*



**Intended use:** NephroPatch™ M is an elastic strip for closing the puncture site after puncturing. The adhesive strip NephroPatch™ M is intended for haemodialysis applications. This external closure strip is used

as a substitute for manual or mechanical compression of the puncture sites by the patient or nursing staff. NephroPatch™ M is a single use aid for use in chronic haemodialysis.